# **Urinalysis Compliance Tools**

POCC Webinar January 19, 2011 Dr. Susan Selgren

### Learning Objectives

- Be able to review and improve upon a laboratory plan for compliance including:
  - Competency
  - Documentation
  - Proficiency
- Understand elements of compliance specific to urinalysis testing
- Gain ability to find resources to help with compliance

#### **Compliance Basics**

- Create your site plan
  - Diagnostic testing
  - Quality control
  - Operator training
- Follow Through!
  - You must actually do what you said you would do!
- Document the implementation
  - To provide evidence that you followed your site's written plan

## Clinical Laboratory Improvement Act (CLIA)

- The Center for Medicare and Medicaid Services (CMS) regulates all laboratory tests in the US
- Regulate through Clinical Laboratory Improvement Amendments (CLIA)
- CLIA Waived labs
  - Perform simple tests
  - Operator might not be a trained laboratory technician
  - □ No legally required audits although audits can occur randomly
- CLIA Moderate Complexity labs
  - Permitted to run more complex tests
  - □ Requires a higher level of operator training / education
  - □ Audits legally required; more frequent

Question for you: CLIA Certification

### **CLIA** Certification

- Waived tests can be performed under a Moderate Complexity certificate
- Moderate Complexity tests cannot be performed under a Waived certificate
- Hospital may hold the CLIA certificate responsible for the testing at offsite facilities
  - outpatient clinics
  - urgent care centers
  - nursing homes
  - physician's office groups



### **CLIA** Waived Tests

- Visual urinalysis testing is waived
- Tests run on instruments can be waived
- Assay menu can include waived and moderately complex tests





### POC Testing Accreditation and Audits

- College of American Pathologists (CAP)
- Commission on Office Laboratory Accreditation (COLA)
- Food and Drug Administration (FDA)
  - CLIA Waived
  - Clinical Laboratory Improvement Act (CLIA) Moderate / High Complexity
    - Moderate Complexity includes Provider Performed Microscopy Procedures (PPMP)
- The Joint Commission accredits non-waived labs; does not recognize waived status
- Survey (audit) every two years to accredit and check compliance

Question for you: Accrediting Agencies

### **Accreditation Tools**



- Example from The Joint Commission is Periodic Performance Review (PPR) software
  - Web-enabled tool via secure extranet
  - Self Assessment Process
  - Submitted annually
  - Plans of Action and Measures of Success
- "Most commonly found deficiencies in laboratory inspections"
  - □ W. Foubister in CAP Today, March 2001 at www.CAP.org

 Each agency has their own auditing process depending upon the situation.



U.S. Food and Drug Administration

- Each agency has their own auditing process depending upon the situation.
- Often the agency leverages an audit checklist



**FD** U.S. Food and Drug Administration

- Each agency has their own auditing process depending upon the situation.
- Often the agency leverages an audit checklist
- Usually the checklist used by your accrediting agency is available to you -- it is helpful for setting up your site's testing process, quality control and operator training procedures



**FD** U.S. Food and Drug Administration

- Each agency has their own auditing process depending upon the situation.
- Often the agency leverages an audit checklist
- Usually the checklist used by your accrediting agency is available to you -- it is helpful for setting up your site's testing process, quality control and operator training procedures
- Tips for survey:
  - □ Know how to access needed records
  - □ Inform surveyor of your staff's schedule
  - Inform staff that they will be asked questions during the survey
  - Encourage open participation



### Now about Urinalysis...

- Standard dip and read urinalysis testing uses strips containing several tests.
- The strips can be read visually against a color chart on the bottle label



 Some urinalysis strip products are good until the expiration date printed on the bottle label and/or carton (so, there is no shorter opened bottle use life)

#### Now about Urinalysis...

- The strips can be read instrumentally, providing more standardization of the reading and documentation of results.
- When urinalysis strips are read instrumentally, the results can be sent into a Hospital or Laboratory Information System (HIS/LIS)



### Urinalysis Results are Semi-Quantitative

Urinalysis strip results are not quantitative

Statistical measures usually not appropriate include:

- □ Linearity
- Precision
- Standard Deviation (SD)
- Coefficient of Variation (CV)

Statistical treatment of semi-quantitative data is unique

### Truth Table Example

	Reference Negative	Reference Positive
Strip Negative	113/125 = 90%	7/150 = 5%
Strip Positive	12/125 = 10%	143/150 = 95%

Leukocytes	Leukocytes Nitrite	Urinary Tract Infection
Nitrite		
Urobilinogen		
Protein		
рН		
Blood		
Specific Gravity		
Ketone		
Bilirubin		
Glucose		

Leukocytes	Leukocytes Nitrite	Urinary Tract Infection
Nitrite	Ηα	Kidney Cancer
Urobilinogen	pH Occult Blood	Kidney Stones Bladder Cancer
Protein		
рН		
Blood		
Specific Gravity		
Ketone		
Bilirubin		
Glucose		





Glucose







### **Urinalysis Results Utility**

#### Urinalysis results are used to detect:

- □ Urinary Tract Infections (UTIs)
- Early chronic kidney disease
- Metabolic dysfunction (Diabetes)
- □ Liver dysfunction

#### Additional popular urine tests include:

- □ hCG
- □ Tablet tests (to detect reducing sugars, bilirubin, and ketone),
- Microscopic analysis
- Visual determination of color and clarity

## Urinalysis Compliance Challenges

- Assuring quality of test results and documentation
- Keeping procedures up to date
- Maintaining competency for large number of operators spread out
- Performing and documenting quality control tests at the required frequencies



## Urinalysis Compliance Challenges

Recording reference ranges in the patient record

- According to The Joint Commission, only quantitative results require this
- Maintaining audit documentation including:
  - Lot numbers
  - Operators
  - Instrument serial numbers
- Knowing what urinalysis testing is being done where, including visual reading
- The operator is often a nurse
  Nurses are not laboratorians





DEPARTMENT OF HEALTH & HUMAN SERVICES



FDA U.S. Food and Drug Administration

**Pre-analytics** 

Sample quality & volume, sample & patient identification

**Pre-analytics** 

Sample quality & volume, sample & patient identification

POC analyzers / Analytics Reagent quality & stability, lock-out features, operatorindependent testing

**Pre-analytics** 

Sample quality & volume, sample & patient identification

POC analyzers / Analytics Reagent quality & stability, lock-out features, operatorindependent testing

Connectivity & Data Management

Result transmission & documentation, remote device management, operator management

**Pre-analytics** 

Sample quality & volume, sample & patient identification

POC analyzers / Analytics Reagent quality & stability, lock-out features, operatorindependent testing

Connectivity & Data Management

Result transmission & documentation, remote device management, operator management

Operator Training & Competency Testing

eLearning, training records, competency testing, autorecertification

**Pre-analytics** 

Sample quality & volume, sample & patient identification

POC analyzers / Analytics Reagent quality & stability, lock-out features, operatorindependent testing

Connectivity & Data Management

Result transmission & documentation, remote device management, operator management

Operator Training & Competency Testing eLearning, training records, competency testing, autorecertification

Inspection Readiness Policies & procedures, safety of testing environment, reagent labeling

### **Quality Control**

- Each laboratory establishes their own quality control process, following manufacturer's instructions.
- Urinalysis strip quality control required for both visual and instrumental testing
- Testing of both normal and abnormal control solutions:
  - Once per day
  - Each time a new lot is used
  - □ And if a new operator is being trained
- Complete quality control records include:
  - Results
  - Date and time
  - Operator
  - Test lot number
  - Quality control lot number
  - Instrument serial number



### Centralized QC management



### **Quality Control Made Easier**

- Instrument reading of urinalysis strips provides:
  - Documentation of quality control results on a printout
  - Quality control results transmission to the HIS / LIS
  - □ Ability to recall QC data for auditor
  - Scheduling of quality control testing
  - Quality control lockout function


### Automated Instrument Quality Checks

- Humidity exposure
- Identify and assure that the correct strip is being used
- Barcode reading of urinalysis strip and quality control material lot numbers
- Automatic calibration on-board
- Combined record of Urinalysis strip and Quality Control lot numbers and expiration dates



# Urinalysis Operators (Users) include:

- Nurses
- Physicians
- Emergency Department (ED) staff
- Physician Assistants
- Nursing Techs
- Nurse Practitioners
- Lab Techs











Operators trained to run each type of diagnostic test







- Operators trained to run each type of diagnostic test
- Operators review training material and take a quiz to show competency







- Operators trained to run each type of diagnostic test
- Operators review training material and take a quiz to show competency
- Supervisor observe testing to assure proper technique







- Operators trained to run each type of diagnostic test
- Operators review training material and take a quiz to show competency
- Supervisor observe testing to assure proper technique
- Documentation of operator certification is key
  - Should include recertification date







- Operators trained to run each type of diagnostic test
- Operators review training material and take a quiz to show competency
- Supervisor observe testing to assure proper technique
- Documentation of operator certification is key
  - Should include recertification date
- Instrument can include an operator lockout function







### Proficiency and Peer Group Testing

- Proficiency testing allows a lab to test unknown solutions and have the results graded
- Proficiency testing programs must be run by a group other than the test manufacturer in order to count for accreditation
- Peer Group testing allows comparison of results from all peer labs
- Programs from diagnostic test manufacturers can help laboratories prepare for accreditation

# **Goals of Competency Assessment**

#### Ensure proficient users

- □ Knowledge and ability to perform a test
- □ Can perform appropriately without supervision
- Determine validity of test results
- □ Take appropriate action when required
- Minimize testing errors
- Ensure regulatory compliance (reduce inspection / survey citations or deficiencies)
- Improve patient safety
- Resource: Some diagnostic test manufacturers offer product education and operator training programs via internet

# How to Organize?

### POCT Coordinators should consider whether to:

- Manage all POC testing under one complexity category (for example, Moderate) to simplify oversight
- Manage POC testing according to the FDA CLIA category that is appropriate

# How to Organize?

### POCT Coordinators should consider whether to:

- Manage all POC testing under one complexity category (for example, Moderate) to simplify oversight
- Manage POC testing according to the FDA CLIA category that is appropriate
- Resource: Clinical Laboratory Standards Institute (CLSI) POCT08-A Quality Practices in Non-Instrumented Near Patient Testing

### Compliance Tips

- Create a program that fits your site's testing and staff resource situation.
- Continuously check on what works, what doesn't.
- Change the program as needed to optimize.
- Keep lines open between Point of Care Coordinators and Nursing Education group.

# Compliance Tips

- Leverage Point of Care testing checklists (CAP and CLIA).
- Watch operators in real time to see who needs training.
- Document competency of staff and audit trail.
- Encourage operators to openly participate in surveys / audits.
- Understand and reduce pre-analytical and data entry error sources.
- Follow manufacturer's instructions for quality control testing and product storage (temperature and humidity).

### Resources

- The Joint Commission
  - □ <u>www.jointcommission.org</u>
- Commission on Office Laboratory Accreditation (COLA)
  - □ <u>www.cola.org</u>
- Centers for Medicare and Medicaid Services (CMS)
  - CLIA: www.cms.hhs.gov/clia/
  - □ CoPs: <u>www.cms.hhs.gov/CFCsandCoPs/</u>
- Centers for Disease Control and Prevention (CDC)
  - www.phppo.cdc.gov/clia/
- Food and Drug Administration CLIA database
  - www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfCLIA/Search.cfm
- College of American Pathology (CAP)
  - □ <u>www.cap.org</u>
- Clinical Laboratory and Standards Institute (CLSI)
  - □ <u>www.clsi.org</u>

### These were today's objectives:

- Be able to review and improve upon a laboratory plan for compliance including:
  - Competency
  - Documentation
  - Proficiency
- Understand elements of compliance specific to urinalysis testing
- Gain ability to find resources to help with compliance